



TIMESCO
ION AIR

ANAESTHETIC & RESUSCITATION
CONSUMABLES

RESUSCITATION

Timesco[™]

ION-AIR RESUSCITATORS - SINGLE PATIENT USE

VERSATILE, RELIABLE & DURABLE

The new Ion-Air manual resuscitators are supplied complete with oxygen tubing, an oxygen reservoir bag and a suitably sized anatomical air cushion mask. The Ion-Air resuscitator is pre-assembled for immediate use and is packed in a tamper proof sealed bag, which is double wrapped in a drawstring bag.

The Ion-Air manual resuscitators are available in a range of sizes, including an infant 300ml, a child 550ml and three adult sizes 1000ml, 1400ml & 1600ml.

The child and infant bags incorporate a 40cm H₂O safety pressure release valve to protect against over inflation of the lungs.

With its textured grip and 360° swivel, the Ion-Air resuscitator is extremely versatile and easy to use with either a face mask or endotracheal tube. The transparent bag offers complete internal visibility.

Key Features:

- Textured grip (non-slip)
- 360° swivel connector for greater versatility
- Glued components (with exception of face mask)
- Adjustable air cushion face mask for maximum patient comfort
- Transparent bag for complete internal visibility
- Excellent recoil recovery following bag compression
- Available in a variety of sizes
- Reliable & durable
- Latex free



NHS CODE	TIMESCO CODE	DESCRIPTION
FDE693	TDM-BVM-1600	ADULT RESUSCITATORS 1600ML SIZE 5 MASK
FDE695	TDM-BVM-1400	ADULT RESUSCITATORS 1400ML SIZE 5 MASK
FDE675	TDM-BVM-1000	ADULT RESUSCITATORS 1000ML SIZE 5 MASK
FDE676	TDM-BVM-0550	CHILD RESUSCITATORS 550ML SIZE 1 & 2 MASKS
FDE694	TDM-BVM-0300	INFANT RESUSCITATORS 300ML SIZE 0 MASK

*Prices are subject to ordering in case quantities.



ION-AIR RESCUER MASKS WITH BACTERIAL & VIRAL FILTER

The Ion-Air Rescuer Mask is a single use emergency mask, which is supplied complete with a one way bacterial and viral filter, built in O₂ port and an elasticated head strap. The mask is supplied in a hard plastic carry case.

NHS CODE	TIMESCO CODE	DESCRIPTION
FDD2552	TDM-PM-250/20	ION-AIR RESCUER MASK IN CASE (PACK OF 20)
	TDM-PM-250	ION-AIR RESCUER MASK IN CASE

TO ORDER VIA NHS SUPPLY CHAIN VISIT WWW.SUPPLYCHAIN.NHS.UK

FACE SHIELD WITH VALVE

KEYRING POUCH

The Timesco single use face shield is designed for convenient, effective protection from cross infection during mouth-to-mouth resuscitation.

The face shield incorporates a one way valve for optimal safety and also features integrated elastic ear bands that fit over the patients ears to keep the mask in place.

Supplied in a handy keyring pouch.



Available in Green or Pink.



NHS CODE	TIMESCO CODE	DESCRIPTION
FDG1065	TDM-PM-050/50	FACE SHIELD IN KEYRING POUCH - GREEN - BOX OF 50
	TDM-PM-055/50	FACE SHIELD IN KEYRING POUCH - PINK - BOX OF 50



FACE SHIELD WITH FILTER

FOIL SACHET

A face shield for personal protection during mouth-to-mouth resuscitation. Supplied in easy tear foil sachet. Single Use.

NHS CODE	TIMESCO CODE	DESCRIPTION
FDG872	TDM-PM-165	FACE SHIELD IN SACHET

RESCUER MVP SUCTION PUMP

LIGHTWEIGHT, HIGH VACUUM PUMP

The Rescuer Manual Vacuum Pump is a lightweight, high vacuum, portable emergency device which can be powered with one hand operation, leaving the other hand free for different tasks. It was designed with a 100° swivel which gives the operator better accessibility whilst working in tight situations and avoids premature engagement of the overflow mechanism. This unit was designed to provide efficient operation and simple maintenance.

The Rescuer pump has multiple settings, allowing the pump to be used from neonates to large adults.



TIMESCO CODE	DESCRIPTION
TDM-MVP-100	RESCUER MVP SUCTION PUMP ASPIRATOR
TDM-MVP-200	SPARE COLLECTION JAR WITH CATHETERS

CF3 FINGER TIP PULSE OXIMETER

FITS A WIDE RANGE OF FINGER SIZES

- 3 year guarantee
- 55mm (w) x 32mm (h) x 28mm (d)
- 30 hours of continuous use
- Large, bright display
- Small & lightweight
- Low battery indicator
- Audible & visual alarm



NHS CODE	TIMESCO CODE	DESCRIPTION
FBC268	TDM-PO-CF3	FINGER TIP PULSE OXIMETER CF3

The CF3 & C52 Fingertip pulse oximeter's are self contained units that are compact and lightweight, designed for spot checking patient SpO2. The units are ideal for use in clinical, ambulance or homecare environments and are suitable for use on both paediatrics and adults.

The displays are OLED, dual colour and provides an easy to read, crisp, bright output that has 4 different direction displays and 6 different display options. Along with showing the usual oxygen levels and pulse rate, the units also display a perfusion indicator, pulse sounds and a Pleth waveform that helps the user identify the quality of the signal that the pulse oximeter is obtaining.

These devices have a low power consumption and can be used continuously for 30 hours. They feature an auto power off design, to conserve battery life and are powered by convenient AAA alkaline batteries with a low battery warning alarm.

C52 FINGER TIP PULSE OXIMETER

DESIGNED WITH PAEDIATRICS IN MIND



Six display modes

- Special design for paediatric & adult use
- Lightweight
- 3 year guarantee
- Display SpO2, PR, pulse bar & plethysmogram
- Low power consumption
- Auto switch off
- Battery low indicator
- Six display modes
- Convenient AAA size alkaline batteries
- Adjustable brightness

NHS CODE	TIMESCO CODE	DESCRIPTION
FBC269	TDM-PO-C52	FINGER TIP PULSE OXIMETER C52

300M HAND HELD PULSE OXIMETER

DEPENDABLE, RUGGED & EASY TO USE

The 300M can continuously monitor, spot check, replay the data in real time and display SpO2 plethysmograph. It can be widely applied in hospitals and clinics.

- Suitable for adult & paediatric use
- High resolution 2.8" TFT display
- Adjustable audio & visual alarms
- Patient ID ranges from 1 to 127
- Auto power off option
- Convenient 2 x AA batteries
- Blue tooth compatible

NHS CODE	TIMESCO CODE	DESCRIPTION
FBC270	TDM-PO-300M	HAND HELD PULSE OXIMETER 300M



SINGLE-USE AIR CUSHION FACE MASKS

MAXIMUM PATIENT COMFORT

Our single-use air cushion face masks are available in a full range of sizes from neonatal to large adult. The air cushion face masks are transparent which allows a clear view of the patient.

Latex free.

NHS CODE	TIMESCO CODE	DESCRIPTION
FDD2570	TDM-MP-1200	AIR CUSHION & HOOK RING, SIZE 0 - BOX OF 50
FDD2569	TDM-MP-1201	AIR CUSHION & HOOK RING, SIZE 1 - BOX OF 50
FDD2565	TDM-MP-1202	AIR CUSHION & HOOK RING, SIZE 2 - BOX OF 50
FDD2566	TDM-MP-1203	AIR CUSHION & HOOK RING, SIZE 3 - BOX OF 32
FDD2567	TDM-MP-1204	AIR CUSHION & HOOK RING, SIZE 4 - BOX OF 32
FDD2568	TDM-MP-1205	AIR CUSHION & HOOK RING, SIZE 5 - BOX OF 32

Our single-use air cushion face masks are also available with an adjustable check valve to ensure minimum leakage and maximum patient comfort.

Latex free.

NHS CODE	TIMESCO CODE	DESCRIPTION
FDD2560	TDM-MP-1210	AIR CUSHION, CHECK VALVE & HOOK RING, SIZE 0 - BOX OF 50
FDD2561	TDM-MP-1211	AIR CUSHION, CHECK VALVE & HOOK RING, SIZE 1 - BOX OF 50
FDD2559	TDM-MP-1212	AIR CUSHION, CHECK VALVE & HOOK RING, SIZE 2 - BOX OF 50
FDD2564	TDM-MP-1213	AIR CUSHION, CHECK VALVE & HOOK RING, SIZE 3 - BOX OF 32
FDD2562	TDM-MP-1214	AIR CUSHION, CHECK VALVE & HOOK RING, SIZE 4 - BOX OF 32
FDD2563	TDM-MP-1215	AIR CUSHION, CHECK VALVE & HOOK RING, SIZE 5 - BOX OF 32



Check Valve

Our single-use air cushion face masks are also available strawberry scented, and with an adjustable check valve to ensure minimum leakage and maximum patient comfort.

Latex free.

NHS CODE	TIMESCO CODE	DESCRIPTION
FDD2556	TDM-MP-1220	AIR CUSHION, CHECK VALVE & HOOK RING SCENTED, SIZE 0 - BOX OF 50
FDD2554	TDM-MP-1221	AIR CUSHION, CHECK VALVE & HOOK RING SCENTED, SIZE 1 - BOX OF 50
FDD2558	TDM-MP-1222	AIR CUSHION, CHECK VALVE & HOOK RING SCENTED, SIZE 2 - BOX OF 50
FDD2555	TDM-MP-1223	AIR CUSHION, CHECK VALVE & HOOK RING SCENTED, SIZE 3 - BOX OF 32
FDD2553	TDM-MP-1224	AIR CUSHION, CHECK VALVE & HOOK RING SCENTED, SIZE 4 - BOX OF 32
FDD2557	TDM-MP-1225	AIR CUSHION, CHECK VALVE & HOOK RING SCENTED, SIZE 5 - BOX OF 32



SILICONE REUSABLE FACE MASKS

FULLY AUTOCLAVEABLE

Silicone face masks are transparent to allow visual checks for bleeding and vomiting. Available in neonate, infant, paediatric, medium, standard and large adult sizes.

Latex free and reusable.



NHS CODE	TIMESCO CODE	DESCRIPTION
FDD668	TDM-MP-1500	SILICONE NEONATAL FACE MASK SIZE 0
FDD669	TDM-MP-1501	SILICONE INFANT FACE MASK SIZE 1
FDD670	TDM-MP-1502	SILICONE PAEDIATRIC FACE MASK SIZE 2
FDD671	TDM-MP-1513	SILICONE SMALL ADULT FACE MASK SIZE 3
FDD1055	TDM-MP-1514	SILICONE MEDIUM ADULT FACE MASK SIZE 4
FDD1056	TDM-MP-1515	SILICONE LARGE ADULT FACE MASK SIZE 5



NASOPHARYNGEAL AIRWAYS

UNIQUE DESIGN

The nasopharyngeal airway has a unique oval flange design, which means the device sits comfortably for the patient, and there is no need for a safety pin. The one piece design is made from a transparent material which enables blockages to be easily and quickly identified.

Latex free. Box of 10.


TIMESCO CODE	DESCRIPTION
TDM-D-106	NASOPHARYNGEAL AIRWAY, SIZE 6 - BOX OF 10
TDM-D-107	NASOPHARYNGEAL AIRWAY, SIZE 7 - BOX OF 10
TDM-D-108	NASOPHARYNGEAL AIRWAY, SIZE 8 - BOX OF 10
TDM-D-109	NASOPHARYNGEAL AIRWAY, SIZE 9 - BOX OF 10

BOUGIE INTRODUCERS

ASSIST DIRECTIONAL CONTROL

The single use bougie introducers are latex free, sterile, individually wrapped and supplied in boxes of 10. Bougie introducers will assist in directional control of the Endotracheal Tube, especially when the laryngeal opening is not fully visible.

Available in neonate, child and adult sizes.



NHS CODE	TIMESCO CODE	DESCRIPTION
FSM742	TDM-D-1130	BOUGIE STRAIGHT 5CH 470MM, NEONATE - BOX OF 10
FSM743	TDM-D-1131	BOUGIE STRAIGHT 10CH 600MM, CHILD - BOX OF 10
	TDM-D-1136	BOUGIE ANGLED 15CH 600MM, ADULT - BOX OF 10
FSM744	TDM-D-1138	BOUGIE ANGLED 15CH 800MM, ADULT - BOX OF 10



Timesco of London Ltd. Timesco House, 3 Carnival Park, Carnival Close, Basildon, Essex SS14 3WN

T: +44 (0) 1268 297 710

F: +44 (0) 1268 297 801

E: sales@timesco.com

www.timesco-primarycare.co.uk

