

COMPATIBILITY CHART



Laryngoscope Blades and Handles

The table below shows what laryngoscope blades are compatible with certain handles. This guide follows the standard light (black) and green light fitting (green) colour coding that are currently used- on all Timesco's laryngoscope systems.

		TIMESCO HANDLE BRAND					
		XENON	XLED	Callisto S	Europa S	Orion	Lunar
TIMESCO BLADE BRAND	Optima CLX	✓	✓	X	X	X	X
	Sirius	✓	✓	X	X	X	X
	Eclipse	✓	✓	X	X	X	X
	Callisto	✓	✓	✓	X	X	X
	Freeway	✓	✓	✓	X	X	X
	Europa	X	X	X	✓	✓	✓

For more information or to order any of the products listed please call 01268 297 700

F.L.022/1